



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse STXBP2 Monoclonal Antibody

单克隆抗体

产品货号: MA02144

第三版

免疫原:Fusion protein of STXBP2

亚型:IgG1

克隆:3322

描述:Syntaxin-binding protein 2 (STXBP2,also known as Unc18-2) belongs to the STXBP/unc-18/SEC1 family. STXBP2 is involved in intracellular trafficking,control of SNARE (soluble NSF attachment protein receptor) complex assembly,and the release of cytotoxic granules by natural killer cells. Mutations in the STXBP2 gene are associated with familial hemophagocytic lymphohistiocytosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF