

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TELO2 Monoclonal Antibody

单克隆抗体

产品货号: MA02148

第三版

免疫原:Fusion protein of TELO2

亚型:IgG1

克隆:950

描述:TELO2 gene encodes a telomere length regulation protein TEL2 homolog. TELO2 may be involved in telomere length regulation and can form TTT complex with TTI1 and TTI2. TTT complex is required to stabilize protein levels of PIKK family proteins and is involved in the cellular resistance to DNA damage stresses,like ionizing radiation (IR),ultraviolet (UV) and mitomycin C (MMC). The activity of mTORC1 and mTORC2 complexes,which regulate cell growth and survival in response to nutrient and hormonal signals,can be promoted,stabalized and maintained by TELO2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF