



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TLE1 Monoclonal Antibody

单克隆抗体

产品货号: MA02152

第三版

免疫原:Fusion protein of TLE1

亚型:IgG2b

克隆:852

描述:Transducin-like enhancer protein 1(TLE1) is a transcriptional corepressor. It can bind to and regress expression of many transcription factors such as NFkB,which is mediated by FOXA2,CTNNB1 and TCF family member in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Expression of TLE1 is highest in skeletal muscle and liver,moderate in placenta,heart,kidney,and spleen,and weak in all other tissues including brain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF