

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TLE3 Monoclonal Antibody

单克隆抗体

产品货号: MA02153

第三版

免疫原:Fusion protein of TLE3

亚型:IgG1

克隆:654

**描述:**TLE3,also named as KIAA1547 and ESG3,belongs to the WD repeat Groucho/TLE family. It is a transcriptional corepressor that binds to a number of transcription factors. It inhibits the transcriptional activation mediated by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. TLE3 is a member of the transducin-like enhancer of split (TLE) family of proteins that have been implicated in the tumorigenesis and classification of sarcomas. It is a candidate biomarker of response to taxane therapy.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF