



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TMEM106B Monoclonal Antibody

单克隆抗体

产品货号: MA02154

第三版

免疫原:Fusion protein of TMEM106B

亚型:IgG1

克隆:123

描述:TMEM106B is a genetic risk factor for frontotemporal lobar degeneration with TDP-43 inclusions (FTLD-TDP). Amyotrophic lateral sclerosis (ALS), like FTLD-TDP, is characterized by pathological TDP-43 inclusions. TMEM106B expression in the brain may be linked to mechanisms of disease in FTLD-TDP and risk alleles confer genetic susceptibility by increasing gene expression.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF