
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TMEM70 Monoclonal Antibody

单克隆抗体

产品货号: MA02155

第三版

免疫原:Fusion protein of TMEM70

亚型:IgG2b

克隆:321

描述:TMEM70 belongs to the TMEM70 family. It is involved in biogenesis of mitochondrial ATP synthase. Defects in TMEM70 are a cause of mitochondrial encephalocardiomyopathy neonatal due to ATP synthase deficiency (MT-ATPSD). A publication identified TMEM70 gene defect as a pan-ethnic disorder and further redefined it as the most common cause of nuclear-origin ATP synthase deficiency.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF