

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human P63 Monoclonal Antibody

单克隆抗体

产品货号: MA02158

第三版

免疫原:Fusion protein of TP63

亚型:IgG1

克隆:447

描述:TP63,also named KET,P63,P73H,P73L and TP73L,belongs to the p53 family. It is a homologue of the tumor suppressor p53 and p73 genes. It is involved in malignancy acquisition and maintenance of cells. Unlike p53,the p63 gene encodes multiple isotypes with remarkably divergent abilities to transactivate p53 reporter genes and induce apoptosis. TP63 acts as a sequence specific DNA binding transcriptional activator or repressor. TP63 has 12 isoforms with MW 40kd(P40),50kd(P60),63kd(P63) and 73kd(P73L). It is Nuclear stain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF