



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TPI1 Monoclonal Antibody

单克隆抗体

产品货号: MA02159

第三版

免疫原:Fusion protein of TPI1

亚型:IgG1

克隆:134

**描述:**TPI1(Triosephosphate isomerase) is also named as TIM,TPI and belongs to the triosephosphate isomerase family. It is a key enzyme in cell metabolism that controls the glycolytic flow and energy production through the interconversion of dihydroxyacetone phosphate (DHAP) and D-glyceraldehyde 3-phosphate (G3P). Defects in TPI1 are the cause of triosephosphate isomerase deficiency (TPI deficiency). It can exist as a homodimer.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF