

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TPP2 Monoclonal Antibody

单克隆抗体

产品货号: MA02160

第三版

免疫原:Fusion protein of TPPII

亚型:IgG2b

克隆:360

描述:TPPII(Tripeptidyl-peptidase 2) is an evolutionary conserved serine peptidase of the subtilisin family that removes tripeptides from the free NH<sub>2</sub> terminus of oligopeptides. The 138-kDa TPP2 polypeptide forms large oligomeric complexes of >1,000 kDa that are detected in the cytoplasm of all cell types. In addition,a membrane associated TPP 2 isoform is present in most cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF