

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TPSAB1 Monoclonal Antibody

单克隆抗体

产品货号: MA02161

第三版

免疫原:Fusion protein of TPSAB1

亚型:IgG1

克隆:790

描述:TPSAB1(Tryptase alpha/beta-1) is also named as TPS1,TPS2,TPSB1 and belongs to the peptidase S1 family. It is the major neutral protease present in mast cells and secreted upon the coupled activation-degranulation response of this cell type. TPSAB1 has 2 isoforms produced by alternative splicing with the molecular weight of 31 kDa and 29 kDa. This full length protein has two glycosylation sites.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF