

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TTF1 Monoclonal Antibody

单克隆抗体

产品货号: MA02163

第三版

免疫原:Synthetic Peptide

亚型:IgG2a

克隆:307

描述:TTF1 is a multifunctional nucleolar protein that terminates ribosomal gene transcription,mediates replication fork arrest and regulates RNA polymerase I transcription on chromatin. It plays a dual role in rDNA regulation,being involved in both activation and silencing of rDNA transcription. Interaction with BAZ2A/TIP5,TTF1 recovers DNA-binding activity. The antibody is specific to TTF1.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF