



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Transferrin Monoclonal Antibody

单克隆抗体

产品货号: MA02164

第三版

免疫原:Fusion protein of Transferrin

亚型:IgG1

克隆:657

描述:Transferrin(TF) is also named as Beta-1 metal-binding globulin,Siderophilin and belongs to the transferrin family. It is a bilobal iron-binding protein that plays the crucial role of binding ferric iron and keeping it in solution,thereby controlling the levels of this important metal. Serum transferrin may have a further role in stimulating cell proliferation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF