



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TAGLN2 Monoclonal Antibody

单克隆抗体

产品货号: MA02165

第三版

免疫原:Fusion protein of Transgelin 2

亚型:IgG1

克隆:370

描述:The transgelin family is a group of proteins that belong to 22kd actin-related corpnin superfamily. Of all three isoforms,transgelin 1 is the best characterized. Transgelin 1,also known as SM22alpha,is a specific marker for differentiated smooth muscle cells. Transgelin 2,also known as SM22 beta,is expressed by both smooth muscle and non-smooth muscle cells in a temporally and spatially regulated pattern. Trangenlin 3,also known as NP25,is only found in highly differentiated neuronal cells. Our test result showed that this antibody is specific to TAGLN2. It does not bind TAGLN1 and TAGLN3.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF