

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat TAGLN Monoclonal Antibody

单克隆抗体

产品货号: MA02166

第三版

免疫原:Fusion protein of transgelin/SM22

亚型:IgG1

克隆:678

描述:The transgelin family is a group of proteins that belong to 22kd actin-related corpnin superfamily. Of all three isoforms,transgelin 1 is the best characterized. Transgelin 1,also known as SM22alpha,is a specific marker for differentiated smooth muscle cells. Transgelin 2,also known as SM22 beta,is expressed by both smooth muscle and non-smooth muscle cells in a temporally and spatially regulated pattern. Trangenlin 3,also known as NP25,is only found in highly differentiated neuronal cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF