

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRPC4AP Monoclonal Antibody

单克隆抗体

产品货号: MA02167

第三版

免疫原:Fusion protein of TRPC4AP

亚型:IgG3

克隆:256

描述:TRPC4AP,also named as C20orf188 and TRRP4AP,is a substrate-specific adapter of a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex required for cell cycle control. TRPC4AP participates in the activation of NFKB1 in response to ligation of TNFRSF1A,possibly by linking TNFRSF1A to the IKK signalosome. The calculated molecular weight of TRPC4AP is 91 kDa. Western blot analysis detected a protein of about 82 kDa expressed in several mouse cell lines and in a human promonocytic cell line.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF