

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human UBE2C Monoclonal Antibody

单克隆抗体

产品货号: MA02172

第三版

免疫原:Fusion protein of UBE2C

亚型:IgG1

克隆:230

描述:UBE2C(Ubiquitin-conjugating enzyme E2 C) is also known as UBCH10 and belongs to the ubiquitin-conjugating enzyme family. It is highly expressed in a variety of human primary tumors,and it is able to promote cell proliferation and malignant transformation. UBE2C,along with the E3 ligase of the anaphasepromoting complex (APC),catalyzes the ubiquitination of mitotic cyclins A and B,as well as securin. UBE2C is essential for cell cycle progression,and mutation of its active site cysteine confers a dominant-negative phenotype.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF