

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human USP13 Monoclonal Antibody

单克隆抗体

产品货号: MA02175

第三版

免疫原:Fusion protein of USP13

亚型:IgG1

克隆:640

描述:USP13(Ubiquitin carboxyl-terminal hydrolase 13) is also named as ISOT3 and belongs to the peptidase C19 family. It is involved in ubiquitin-dependent protein catabolism. USP13 is responsible for MITF deubiquitination,utilizing a short hairpin RNA library against known deubiquitinating enzymes. This protein has 2 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF