

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human VAPB Monoclonal Antibody

单克隆抗体

产品货号: MA02176

第三版

免疫原:Fusion protein of VAPB

亚型:IgG1

克隆:437

描述: Vesicle-associated membrane protein-associated protein B (VAPB) is an integral membrane protein localized to the endoplasmic reticulum (ER) membrane. VAPB has been implicated in various cellular processes, including ER stress, the unfolded protein response (UPR) and calcium homeostasis regulation. The mutations in the gene of VAPB cause amyotrophic lateral sclerosis 8 (ALS8) and some other related forms of motor neuron disease including late onset spinal muscular atrophy.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF