



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat 14-3-3 theta Monoclonal Antibody

单克隆抗体

产品货号: MA02186

第三版

免疫原:Fusion protein of 14-3-3

亚型:IgG1

克隆:954

描述:14-3-3 proteins interact with a wide spectrum of proteins and possess diverse functions. Mammals express seven distinct 14-3-3 isoforms (gamma,epsilon,beta,zeta,sigma,theta,tau) that form multiple homo- and hetero- dimmers. 14-3-3 proteins display the highest expression levels in the brain,and have been implicated in several neurodegenerative diseases,including Alzheimer's disease and amyotrophic lateral sclerosis. This antibody recognizes 14-3-3 tau (YWHAQ),eta (YWHAH),and zeta (YWHAZ),but not epsilon (YWHAE).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF