

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AIRE Monoclonal Antibody

单克隆抗体

产品货号: MA02188

第三版

免疫原:Fusion protein of AIRE

亚型:IgG1

克隆:954

描述:This gene encodes a transcriptional regulator that forms nuclear bodies and interacts with the transcriptional coactivator CREB binding protein. The encoded protein plays an important role in immunity by regulating the expression of autoantigens and negative selection of autoreactive T-cells in the thymus. Mutations in this gene cause the rare autosomal-recessive systemic autoimmune disease termed autoimmune polyendocrinopathy with candidiasis and ectodermal dystrophy (APECED).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF