



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AKAP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02189

第三版

免疫原:Fusion protein of AKAP3

亚型:IgG2a

克隆:79

描述:AKAP3 is a member of A-kinase anchoring proteins (AKAPs), a family of functionally related proteins that target protein kinase A to discrete locations within the cell. AKAP3 is reported to participate in protein-protein interactions with the R-subunit of the protein kinase A as well as sperm-associated proteins. AKAP3 is expressed in spermatozoa and localized to the acrosomal region of the sperm head as well as the length of the principal piece. It may function as a regulator of motility, capacitation, and the acrosome reaction.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF