



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANXA4 Monoclonal Antibody

单克隆抗体

产品货号: MA02190

第三版

免疫原:Fusion protein of Annexin IV

亚型:IgG2c

克隆:3476

描述:Annexin IV /placental anticoagulant protein II (ANX4) is a member of the annexin or lipocortin family of calcium-dependent phospholipid binding proteins, and is located at 2p13. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. ANX4 shares 45 to 59% identity with other members and possibly interacts with ATP. ANX4 has in vitro anticoagulant activity and also inhibits phospholipase A2 activity. ANX4 is almost exclusively expressed in epithelial cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF