

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANXA5 Monoclonal Antibody

单克隆抗体

产品货号: MA02191

第三版

免疫原:Fusion protein of Annexin V

亚型:IgG1

克隆:9540

描述:Annexin A5 (ANXA5), is a member of the structurally related family of annexin proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. Annexin 5 has also been described as placental anticoagulant protein I, vascular anticoagulant-alpha, endonexin II, lipocortin V, placental protein 4 and anchorin CII.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF