

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOA4 Monoclonal Antibody

单克隆抗体

产品货号: MA02192

第三版

免疫原:Fusion protein of APOA4

亚型:IgG1

克隆:4601

描述:APOA4 is a major protein component of chylomicrons in postprandial lymph and plasma. It is synthesized in the small intestine and the liver, attached to chylomicrons by enterocytes, and secreted into intestinal lymph during fat absorption. APOA4 is intimately involved in metabolism and has protective roles in cardiovascular diseases like atherosclerosis. APOA4 can promote insulin secretion and maintain glucose homeostasis. It also has potent antioxidant and anti-inflammatory properties.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF