



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat BDNF Monoclonal Antibody

单克隆抗体

产品货号: MA02194

第三版

免疫原:Fusion protein of BDNF

亚型:IgG2a

克隆:1201

描述:BDNF is a member of the nerve growth factor family. It is induced by cortical neurons, and is necessary for survival of striatal neurons in the brain. Expression of BDNF is reduced in both Alzheimer's and Huntington disease patients. It plays a role in the regulation of stress response and in the biology of mood disorders.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF