



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PAX5 Monoclonal Antibody

单克隆抗体

产品货号: MA02198

第三版

免疫原:Fusion protein of BSAP/PAX5

亚型:IgG1

克隆:96

描述:PAX5,also named B-cell specific activator protein (BSAP),is a member of the highly conserved paired box (PAX)-domain family of transcription factors. PAX5 is essential in normal B-cell lymphopoiesis. It is expressed in pro-B,pre-B and mature B lymphocytes but not in plasma cells. PAX5 is a common target in leukemogenesis of B-ALL (B-lineage acute lymphoblastic leukemia) along with CDKN2A. Owing to its frequent deletion in B-ALL,PAX5 could be used as one of the molecular markers in diagnosis and monitoring of the disease. This antibody is specific to PAX5. It will not cross react with other PAX members.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF