

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CA9 Monoclonal Antibody

单克隆抗体

产品货号: MA02199

第三版

免疫原:Fusion protein of CA9

亚型:IgG2a

克隆:326

描述:CA9(Carbonic anhydrase 9) may be involved in the control of cell proliferation and transformation and appears to be a novel specific biomarker for a cervical neoplasia.It is a tumor-associated antigen that has been shown to have diagnostic utility in identifying cervical dysplasia and carcinoma.The protein is presentboth on the plasma membrane and in the nucleus of cells and has the molecular.It was observed the 54KDa and 58 Kda bands in A549(JOC,vol21,no.3,473-482).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF