



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD200 Monoclonal Antibody

单克隆抗体

产品货号: MA02202

第三版

免疫原:Fusion protein of CD200

亚型:IgG1

克隆:657

描述:This gene encodes a type I membrane glycoprotein containing two extracellular immunoglobulin domains, a transmembrane and a cytoplasmic domain. This gene is expressed by various cell types, including B cells, a subset of T cells, thymocytes, endothelial cells, and neurons. The encoded protein plays an important role in immunosuppression and regulation of anti-tumor activity. Alternative splicing results in multiple transcript variants encoding different isoforms.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF