



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human P-Cadherin Monoclonal Antibody

单克隆抗体

产品货号: MA02205

第三版

免疫原:Fusion protein of CDH3

亚型:IgG1

克隆:9987

描述: Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. P-cadherin (placental cadherin), also known as CDH3 (cadherin 3), is a classical member of the cadherin superfamily which also include E-, N-, R-, and B-cadherins. Aberrant expression of P-cadherin has been reported in many cancers, including breast, gastric, endometrial, ovarian, and prostate cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF