



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CDK6 Monoclonal Antibody

单克隆抗体

产品货号: MA02206

第三版

免疫原:Fusion protein of CDK6

亚型:IgG1

克隆:32

描述:CDK6(Cyclin-dependent kinase 6) is involved in initiation and maintenance of cell cycle exit during cell differentiation and it prevents cell proliferation and regulates negatively cell differentiation,but is required for the proliferation of specific cell types.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF