



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CGB3 Monoclonal Antibody

单克隆抗体

产品货号: MA02207

第三版

免疫原:Fusion protein of CGB/hCG

亚型:IgG2a

克隆:1102

描述:CGB is a member of the glycoprotein hormone beta chain family and encodes the beta 3 subunit of chorionic gonadotropin (CG). Glycoprotein hormones are heterodimers consisting of a common alpha subunit and an unique beta subunit which confers biological specificity. CG is produced by the trophoblastic cells of the placenta and stimulates the ovaries to synthesize the steroids that are essential for the maintenance of pregnancy.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF