

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human CTAGE1 Monoclonal Antibody

单克隆抗体

产品货号: MA02212

---

第三版

免疫原:Fusion protein of CTAGE1

亚型:IgG2b

克隆:55

描述:CTAGE1 (Cutaneous T-Cell Lymphoma-Associated Antigen 1) is a Protein Coding gene. Diseases associated with CTAGE1 include Cutaneous T Cell Lymphoma and Lymphoma. An important paralog of this gene is CTAGE5.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF