



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DAO Monoclonal Antibody

单克隆抗体

产品货号: MA02215

第三版

免疫原:Fusion protein of DAO

亚型:IgG1

克隆:367

描述:D-amino acid oxidase (DAO),with MW of 39 kDa,is a flavoprotein which uses flavin adenine dinucleotide as its prosthetic group. Its substrates include a wide variety of D-amino acids,but it is inactive on the naturally occurring L-amino acids. DAO may play a role in the pathophysiology of schizophrenia and act as a detoxifying agent which removes D-amino acids that accumulate during aging.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF