



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse BRINP1 Monoclonal Antibody

单克隆抗体

产品货号: MA02216

第三版

免疫原:Fusion protein of DBC1

亚型:IgG2c

克隆:206

描述:DBC1,also named as DBCCR1 and FAM5A,inhibits cell proliferation by negative regulation of the G1/S transition. It mediates cell death which is not of the classical apoptotic type and regulates expression of components of the plasminogen pathway.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF