

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ECHS1 Monoclonal Antibody

单克隆抗体

产品货号: MA02217

第三版

免疫原:Fusion protein of ECHS1

亚型:IgG2b

克隆:666

描述:Enoyl-coenzyme A hydratase (ECHS1) is a mitochondrial protein which catalyzes the hydration of 2-trans-enoyl-coenzyme A intermediates to L-3-hydroxyacyl-coenzyme A, playing key role in metabolism of fatty acids in mitochondria. ECHS1 is highly expressed in muscle, liver and fibroblasts. Altered expression of ECHS1 has been considered as a characteristic feature of mitochondria dysfunction. (23275097,23235493)

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF