

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EEA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02218

第三版

免疫原:Synthetic Peptide

亚型:IgM

克隆:777

描述:Early endosome antigen 1 (EEA1) is a marker of early endosomes and has an important role in endosomal trafficking. EEA1 contains coiled-coil regions and an FYVE-type zinc finger which interacts with phosphatidylinositol-3-phosphate. EEA1 functions as a Rab5 effector,mediates endosome docking and,together with SNAREs,leads to membrane fusion.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF