

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ERCC1 Monoclonal Antibody

单克隆抗体

产品货号: MA02220

第三版

免疫原:Fusion protein of ERCC1

亚型:IgG1

克隆:654

描述:Excision Repair Cross Complementing 1 (ERCC1) is a structure-specific endonuclease that is responsible for the 5'-incision during DNA repair. It forms a complex with ERCC11,XPF and ERCC4,which are required in both recombinatorial repair and nucleotide excision repair. It has been found that ERCC1,together with RRM1,are determinants of survival after surgical treatment of early-stage,non-small-cell lung cancer. (Ref: Simon,ER. 2005)

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF