

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FABP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02221

第三版

免疫原:Fusion protein of FABP3

亚型:IgG2a

克隆:156

描述:FABP3 (fatty-acid-binding protein 3),also known as heart-type FABP or mammary-derived growth inhibitor (MDGI),is a small 15-kDa cytoplasmic protein transporting fatty acids and other lipophilic substances from the cytoplasm to the nucleus. It is most ubiquitously expressed in heart and skeletal muscle.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF