



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat FAK Monoclonal Antibody

单克隆抗体

产品货号: MA02223

第三版

免疫原:Fusion protein of FAK

亚型:IgG2c

克隆:304

描述:FAK (Focal adhesion kinase 1) is also named as FAK1,FADK,pp125FAK,FAK and belongs to the protein kinase superfamily. It is a critical tyrosine kinase that modulates cell adhesion,migration,proliferation and survival in response to extracellular signals . It also acts as a pivotal signal 'integrator',controlling and coordinating cellular responses that include cell migration,survival,proliferation and,epithelial tissue repair after DNA damage . This protein has 6 isoforms produced by alternative promoter usage and alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF