

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FUT8 Monoclonal Antibody

单克隆抗体

产品货号: MA02225

第三版

免疫原:Fusion protein of FUT8

亚型:IgG2a

克隆:980

描述:Fucosyltransferase 8 (FUT8) is GDP-fucose-glycoprotein fucosyltransferase and is localized in Golgi apparatus. FUT8 is involved in protein glycosylation by catalyzing the addition of fucose in alpha 1-6 linkage to the first GlcNAc residue, next to the peptide chains in N-glycans. FUT8 may be involved in human aging by regulating protein N-glycosylation and in the molecular pathogenesis of human pulmonary emphysema (PE).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF