

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GOT1 Monoclonal Antibody

单克隆抗体

产品货号: MA02227

第三版

免疫原:Fusion protein of GOT1

亚型:IgG1

克隆:620

描述:Glutamate oxaloacetate transaminase 1 (GOT1) catalyzes the reversible reaction of L-aspartate and alpha-ketoglutarate into oxaloacetate and L-glutamate and plays a key role in carbon and nitrogen metabolism. GOT1 can potentially control the intracellular levels of reactive oxygen species (ROS) through NADPH synthesis and enhances tumor growth. GOT1 expression correlates with the growth of several tumors.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF