

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GSTO1 Monoclonal Antibody

单克隆抗体

产品货号: MA02228

第三版

免疫原:Fusion protein of GSTO1

亚型:IgG1

克隆:165

描述:Glutathione-S-transferase omega 1 (GSTO1) belongs to a new subfamily of GSTs and is the rate limiting enzyme for biotransformation of inorganic arsenic,environmental carcinogen. Expression of GSTO1-1 was abundant in a wide range of normal tissues,particularly liver,macrophages,glial cells,and endocrine cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF