



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HEXIM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02231

第三版

免疫原:Fusion protein of HEXIM1

亚型:IgM

克隆:316

描述:HEXIM1,also named CLP1,EDG1,HIS1 or MAQ1,is best known as the inhibitor of positive transcription elongation factor b,which regulates the transcription elongation of RNA polymerase II and controls 60-70% of mRNA synthesis. HEXIM1 regulates NF-kappa-B,ESR1,NR3C1 and CIITA-dependent transcriptional activity. Previous studies showed that HEXIM1 interacts with two key p53 regulators,nucleophosmin and HDM2,implying a possible connection between HEXIM1 and the p53 signaling pathway.The mainstream MW of this protein is 70 and 54 kDa. This antibody recognises 54 kDa band.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF