



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL1F10 Monoclonal Antibody

单克隆抗体

产品货号: MA02235

第三版

免疫原:Fusion protein of IL1F10

亚型:IgG1

克隆:345

描述:IL1F10 is a member of the interleukin 1 cytokine family. This cytokine is thought to participate in a network of interleukin 1 family members to regulate adapted and innate immune responses.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF