

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL1F5 Monoclonal Antibody

单克隆抗体

产品货号: MA02236

第三版

免疫原:Fusion protein of IL1F5

亚型:IgG2b

克隆:6663

描述:IL36RN is a member of the interleukin 1 cytokine family. This cytokine was shown to specifically inhibit the activation of NF-kappaB induced by interleukin 1 family,member 6 (IL1F6). It inhibits the activity of interleukin-36 (IL36A,IL36B and IL36G) by binding to receptor IL1RL2 and preventing its association with the coreceptor IL1RAP for signaling.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF