



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL36B Monoclonal Antibody

单克隆抗体

产品货号: MA02237

第三版

免疫原:Fusion protein of IL1F8

亚型:IgG1

克隆:456

描述:The protein encoded by this gene is a member of the interleukin 1 cytokine family. Protein structure modeling indicated that this cytokine may contain a 12-stranded beta-trefoil structure that is conserved between IL1A (IL-A alpha) and IL1B (IL-1 beta). This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. Two alternatively spliced transcript variants encoding distinct isoforms have been reported. IL36B (Interleukin 36, Beta) is a Protein Coding gene. Among its related pathways are Innate Immune System and Toll-Like receptor Signaling Pathways. GO annotations related to this gene include cytokine activity and interleukin-1 receptor binding. An important paralog of this gene is IL36A.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF