

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IMP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02239

第三版

免疫原:Fusion protein of IMP3

亚型:IgG2b

克隆:666

描述:This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain. IMP3 (IMP3, U3 Small Nucleolar Ribonucleoprotein) is a Protein Coding gene. Diseases associated with IMP3 include Enchondroma. Among its related pathways are MAPK-Erk Pathway and rRNA processing in the nucleus and cytosol. GO annotations related to this gene include poly(A) RNA binding and structural constituent of ribosome.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF