



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IVD Monoclonal Antibody

单克隆抗体

产品货号: MA02241

第三版

免疫原:Fusion protein of IVD

亚型:IgG2a

克隆:999

描述:IVD(Isovaleryl CoA dehydrogenase,mitochondrial) is a member of the acylCoA dehydrogenase family and is involved in the catabolism of leucine. In size,IVD precursor and mature proteins produced by class I mutants are indistinguishable from their normal counterparts. Class II,III,and IV mutants make IVD precursor proteins 42 kD in size rather than the normal 45 kD. It has 2 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF