

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PSA Monoclonal Antibody

单克隆抗体

产品货号: MA02242

第三版

免疫原:Fusion protein of KLK3/PSA

亚型:IgG1

克隆:456

描述:Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. kallikrein-related peptidase 3, which is also called PSA (prostate-specific antigen) in the clinical setting, is useful in the diagnosis and monitoring of prostatic carcinoma. KLK3 has 5 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF