



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human L2HGDH Monoclonal Antibody

单克隆抗体

产品货号: MA02243

第三版

免疫原:Fusion protein of L2HGDH

亚型:IgG2a

克隆:

描述:2HGDH(L-2-hydroxyglutarate dehydrogenase,mitochondrial) is also named as duranin,C14orf160 and belongs to the L2HGDH family. The putative L2HGDH is predicted to be targeted to the mitochondria where its mitochondrial targeting sequence is presumably removed. Defects in L2HGDH are the cause of L-2-hydroxyglutaric aciduria (L2HGA).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF